WEGO賦高


- Personalized design
- Safe and reliable
- NMPA certification


Single-use Sterile Oviposition Needles

ULTIMATE QUALITY WEGO INTELLIGENT MANUFACTURING

LEADING THE INDUSTRY STANDARD, RADIATION STERILIZATION WITH NO RESIDUE IS SAFER

The first drafting unit of the industry standard YY/T 1698-2020 Medical devices for assisted reproductive technology in vitro - Assisted reproduction aspiration needles





Single chamber ovum retrieval needle

Double chamber ovum retrieval needle



The pen-type needle point is sharp, Reduced tissue damage



Innovative MARK markers, Clear development under B-ultrasound



Special treatment for inner and outer wall of steel needle to reduce surface roughness, inner wall Ra≤0.17 and outer wall Ra≤0.09



Change the end of the suction tube to a straight tube, which can be used for recoil puncture



All catheters use super-slip material to reduce friction damage



The handle is integrally molded to guarantee air tightness and stabilize the pressure in the lumen

Product parameters

Model	Gauge	MARK Area Shape And Length	MARK Area And Needle Tip Distance	Puncture Needle Length (cm)	Suction Catheter Length (cm)	Vacuum Catheter Length (cm)	Flush Catheter Length (cm)	Handle Color
	16 16 16	Punctate 10mm	3mm	35	65 80 100	45	_	Blue
Type I (Single Chamber)	17 17 17 17	Punctate 10mm	3mm	35	55 65 80 100	45	-	Red Red
	18 18 18	Punctate 10mm	3mm	35	65 80 100	45	_	Pink Pink
	16 16 16	Threaded shape 10mm	3mm	35	65 80 100	45	_	Blue
Type III (Single Chamber)	17 17 17	Threaded shape 10mm	3mm	35	55 65 80 100	45	-	Red
	18 18 18	Threaded shape 10mm	3mm	35	65 80 100	45	-	Pink
Type II (Double Chamber)	17 17 17 17	Threaded shape 10mm	3mm	35	55 65 80 100	45	55 65 80 100	Red





Shandong Wego Ruisheng Medical Device Co., Ltd

Address: No. 1, Weigao Road, Torch High-tech Industrial Development Zone, Weihai City, Shandong Province

Weihai City, Shandong Province
Telephone :400-646-6666 Postal code : 264210

Tel: 0631-5715912 Website: www.wegoruisheng.com Fax: 0631-5716403 Email: rsxshscb@weigaogroup.com